AgroEcoLab @ UMD | <u>www.agroecologylab.com</u> Last updated July 3, 2018

Gravimetric Soil Moisture Procedure

Materials:

- Aluminum weigh boats (tins)
- Analytical balance
- Drying oven set to 105 °C
- Sample spoon
- Pencil
- Lab Notebook

Procedure:

- 1. Place a labeled tin on the balance.
- 2. Record the tin number and the tin weight.
- 3. Add 50g of soil to the tin and record the exact weight. You only need to do this once per sample.
- 4. Place the soils in the oven at 105 °C for 48 hours. Note: You can let the soils sit for a few days as long as you recorded the wet weight, first there is no rush to get them in the oven.
- 5. Once the soils are dry, re-weigh the tins plus dry soil. Record the weight.

Calculations and example spreadsheet:

ID	Tin Weight (g)	Wet Weight (g)	Wet Weight + Tin (g)	Dry Weight Soil + Tin (g)	Dry Weight (g)	Grav. %wet wt
001	5.351	50.363	55.714	43.694	38.343	23.9
002	5.351	50.416	55.767	43.794	38.443	23.7
003	5.351	50.608	55.959	43.982	38.631	23.7

1. Dry wt = (tin wt + dry soil wt) - tin wt

2. Gravimetic soil moisture (%) = [(wet wt - dry wt)/ wet wt] *100