

Protocol for measuring cation exchange capacity (CEC) in soils with high calcium content, such as limed, agricultural soils

Purpose

To measure exchangeable cations and cation exchange capacity.

Reagents

- 0.1M barium chloride BaCl_2 (1L)
 - 20.82 g of BaCl_2 in 1L DI water
- 0.002M barium chloride BaCl_2 (1L)
 - 0.416g of BaCl_2 in 1L DI water
- 0.005M* magnesium sulfate heptahydrate $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ (1L)
 - 0.602g $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ in 1L DI water

*A 0.01M solution instead (1.2g $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ /1L DI water) also works, just adjust the calibration curve

Procedure

- Write tube number on data sheet and record weight of tube. Tare balance.
- Weigh 2.00g (± 0.01 g) of soil in empty centrifuge tube. Record weight.
- Add 20.0 ml of 0.1M BaCl_2 .
- Cap tube. Shake for 1 hour.
- Centrifuge tubes at 4000 x g (about 6000 rpm on VWR benchtop centrifuges) for 30 minutes (remember to balance the rotor).
- Decant supernatant liquid (may be saved for exchangeable cation analysis/soil pH measurement). You may use a transfer pipette to remove the last few milliliters of supernatant if needed.
- Add 20.0mL of 0.002M BaCl_2 . Disperse solids with vortex mixer. Shake for 10 minutes.
- Centrifuge for 30 minutes. Decant and discard supernatant.
- Repeat steps 6, 7, and 8 two times.
- After decanting 0.002M BaCl_2 for the third time, weigh the tube plus soil plus entrained solution.
- Add 10.0mL 0.005M $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$. Disperse solids. Shake 30 minutes. NOTE: For critical determinations, the solution ionic strength can be adjusted at this point by matching its electrical conductivity to that of a reference solution. See Gillman 1979 or Rhoades 1982 for a description.
- Let tubes stand overnight, with occasional hand shaking. This long period of equilibration is necessary to allow the Mg^{2+} ions to completely replace the Ba^{2+} ions on the exchange complex. Displaced barium ions precipitate as BaSO_4 .
- Centrifuge 10 minutes. Filter supernatant into a labeled bottle or scintillation vial using Whatman #42 filters.
- Store in the freezer until ready for analysis.

Analysis

Analyze on atomic absorption spectrophotometer for Mg using the Air-Acetylene flame. Standards can be made in DIW and with the 500 mg Mg/L stock solution with the following standard curve. Run a reagent blank as DIW. Subtract out the concentration of $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ in solution blank from the sample concentrations.

Calculations

The total CEC is equivalent to the mass (in millequivalents) of Mg^{++} removed from solution divided by the soil weight.

1. Calculate entrained solution volume:
Entrained = (weight of tube + soil + entrained) - (weight of empty tube + **oven-dry** soil). This number represents the **mass**; assuming a specific gravity of 1 (close enough), it is also the **volume**, in mL.
2. Convert mg/L to meq Mg^{++} in the extract, for both the blank and the sample (remember, $\text{Mass} = \text{Concentration} \times \text{Volume}$; this step also incorporated unit conversions for Mg meq and mL to L).

$$\text{meq Mg}^{++} = \frac{\text{mg Mg}}{\text{L}} \times \frac{1 \text{mmol Mg}}{24.305 \text{ mg Mg}} \times \frac{2 \text{ meq Mg}^{++}}{1 \text{mmol Mg}} \times (10 \text{mL} + \text{entrained}) \times \frac{1 \text{L}}{1000 \text{mL}}$$

3. The difference between meq Mg^{++} in the blank and in the sample represents the mass of magnesium ions adsorbed onto the soil exchange complex. Divide by (dry) soil mass, convert meq to cmol (+) and g to kg, and you have CEC.

CEC (cmol(+)/kg soil)=

$$\frac{(\text{meq Mg}^{++}_{\text{blank}} - \text{meq Mg}^{++}_{\text{sample}})}{2 \text{g moist soil}} \times \frac{(1 + w_d) \text{g moist soil}}{\text{g dry soil}} \times \frac{1 \text{ cmol (+)}}{10 \text{ meq Mg}^{++}} \times \frac{1000 \text{ g soil}}{\text{kg soil}}$$

References

- Gillman, G.P. 1979. A proposed method for the measurement of exchange properties of highly weathered soils. *Aust. J. Soil Res.* 17:129-139.
- Hendershot, W.H., and M. Duquette. 1986. A simple barium chloride method for determining cation exchange capacity and exchangeable cations. *Soil Sci. Soc. Am. J.* 50:605-608.
- Rhoades, J.D. 1982. Soluble Salts. In A.L. Page et al. (ed.) *Methods of soil analysis. Part 2. Agronomy* 9: 167-178.